



ACRO BIOTECH, Inc.

PRODUCT CATALOGUE



CONTENT

Drugs of Abuse 1-4

Cardiac Markers 5

Tumor Markers 5

Women's Health 6-7

Infectious Disease 7-8

Others 9

Specimen: Urine**Acetaminophen (ACE)**

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DAC-102	Urine	Cassette	25 T	5000 ng/mL*	CE
DAC-114	Urine	Panel	25 T		CE
DAC-101	Urine	Dipstick	50 T		CE

7-Aminoclonazepam (ACL)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DACL-102	Urine	Cassette	25 T	300 ng/mL*	CE
DACL-114	Urine	Panel	25 T		CE
DACL-101	Urine	Dipstick	50 T		CE

Alcohol

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DAL-102	Urine	Cassette	25 T	0.02%	CE
DAL-114	Urine	Panel	25 T		CE
DAL-101	Urine	Dipstick	50 T		CE

Amphetamine (AMP) FDA 510(k) Clearances

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DAM-102	Urine	Cassette	25 T	300 ng/mL, 500 ng/mL, 1000 ng/mL*	CE
DAM-114	Urine	Panel	25 T		CE
DAM-101	Urine	Dipstick	50 T		CE

Barbiturate (BAR)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DBA-102	Urine	Cassette	25 T	200 ng/mL, 300 ng/mL*	CE
DBA-114	Urine	Panel	25 T		CE
DBA-101	Urine	Dipstick	50 T		CE

Buprenorphine (BUP)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DBU-102	Urine	Cassette	25 T	5 ng/mL, 10 ng/mL*	CE
DBU-114	Urine	Panel	25 T		CE
DBU-101	Urine	Dipstick	50 T		CE

Benzodiazepines (BZO) FDA 510(k) Clearances

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DBZ-102	Urine	Cassette	25 T	100 ng/mL, 200 ng/mL, 300 ng/mL*, 500 ng/mL	CE
DBZ-114	Urine	Panel	25 T		CE
DBZ-101	Urine	Dipstick	50 T		CE

Clonazepam (CLO)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DCL-102	Urine	Cassette	25 T	150 ng/mL, 400 ng/mL*	CE
DCL-114	Urine	Panel	25 T		CE
DCL-101	Urine	Dipstick	50 T		CE

Cocaine (COC) FDA 510(k) Clearances

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DCO-102	Urine	Cassette	25 T	100 ng/mL, 150 ng/mL, 300 ng/mL*	CE
DCO-114	Urine	Panel	25 T		CE
DCO-101	Urine	Dipstick	50 T		CE

Cotinine (COT)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DCT-102	Urine	Cassette	25 T	10 ng/mL, 50 ng/mL, 100 ng/mL, 200 ng/mL*	CE
DCT-114	Urine	Panel	25 T		CE
DCT-101	Urine	Dipstick	50 T		CE

Diazepam (DIA)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DDI-102	Urine	Cassette	25 T	300 ng/mL*	CE
DDI-114	Urine	Panel	25 T		CE
DDI-101	Urine	Dipstick	50 T		CE

Ethylenediamine-dimethylphosphinic acid (EDDP)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DED-102	Urine	Cassette	25 T	100 ng/mL, 300 ng/mL*	CE
DED-114	Urine	Panel	25 T		CE
DED-101	Urine	Dipstick	50 T		CE

* Standard Cut-off



Drugs of Abuse

Ethyl Glucuronide (ETG)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DET-102	Urine	Cassette	25 T	500 ng/mL*, 1000 ng/mL	CE
DET-114	Urine	Panel	25 T		CE
DET-101	Urine	Dipstick	50 T		CE

Fentanyl (FYL)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DFY-102	Urine	Cassette	25 T	10 ng/mL, 20 ng/mL*	CE
DFY-114	Urine	Panel	25 T		CE
DFY-101	Urine	Dipstick	50 T		CE

Ketamine (KET)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DKE-102	Urine	Cassette	25 T	100 ng/mL*, 300 ng/mL, 500 ng/mL, 1000 ng/mL	CE
DKE-114	Urine	Panel	25 T		CE
DKE-101	Urine	Dipstick	50 T		CE

Lysergic Acid Diethylamide (LSD)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DLS-102	Urine	Cassette	25 T	20 ng/mL*	CE
DLS-114	Urine	Panel	25 T		CE
DLS-101	Urine	Dipstick	50 T		CE

6-Monoacetylmorphine (6-MAM)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DMA-102	Urine	Cassette	25 T	10 ng/mL*	CE
DMA-114	Urine	Panel	25 T		CE
DMA-101	Urine	Dipstick	50 T		CE

Tenamfetamine (MDA)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DMDA-102	Urine	Cassette	25 T	500 ng/mL*	CE
DMDA-114	Urine	Panel	25 T		CE
DMDA-101	Urine	Dipstick	50 T		CE

Ecstasy (MDMA) **FDA 510(k) Clearances**

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DMD-102	Urine	Cassette	25 T	500 ng/mL*, 1000 ng/mL	CE
DMD-114	Urine	Panel	25 T		CE
DMD-101	Urine	Dipstick	50 T		CE

Methylenedioxypyrovalerone (MDPV)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DMDP-102	Urine	Cassette	25 T	1000 ng/mL*	CE
DMDP-114	Urine	Panel	25 T		CE
DMDP-101	Urine	Dipstick	50 T		CE

Methamphetamine (MET) **FDA 510(k) Clearances**

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DME-102	Urine	Cassette	25 T	300 ng/mL, 500 ng/mL, 1000 ng/mL*	CE
DME-114	Urine	Panel	25 T		CE
DME-101	Urine	Dipstick	50 T		CE

Morphine (MOP)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DMO-102	Urine	Cassette	25 T	100 ng/mL, 300 ng/mL*	CE
DMO-114	Urine	Panel	25 T		CE
DMO-101	Urine	Dipstick	50 T		CE

Methylphenidate(MPD)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DMP-102	Urine	Cassette	25 T	1000 ng/mL*, 300 ng/mL	CE
DMP-114	Urine	Panel	25 T		CE
DMP-101	Urine	Dipstick	50 T		CE

Methaqualone (MQL)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DMQ-102	Urine	Cassette	25 T	300 ng/mL*	CE
DMQ-114	Urine	Panel	25 T		CE
DMQ-101	Urine	Dipstick	50 T		CE

Methadone (MTD)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DMT-102	Urine	Cassette	25 T	200 ng/mL, 300 ng/mL*	CE
DMT-114	Urine	Panel	25 T		CE
DMT-101	Urine	Dipstick	50 T		CE

* Standard Cut-off



Opiates (OPI) FDA 510(k) Clearances

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DOP-102	Urine	Cassette	25 T	2000 ng/mL*	CE
DOP-114	Urine	Panel	25 T		CE
DOP-101	Urine	Dipstick	50 T		CE

Oxycodone (OXY)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DOX-102	Urine	Cassette	25 T	100 ng/mL*	CE
DOX-114	Urine	Panel	25 T		CE
DOX-101	Urine	Dipstick	50 T		CE

Phencyclidine (PCP) FDA 510(k) Clearances

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DPC-102	Urine	Cassette	25 T	25 ng/mL*	CE
DPC-114	Urine	Panel	25 T		CE
DPC-101	Urine	Dipstick	50 T		CE

Propoxyphene (PPX)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DPP-102	Urine	Cassette	25 T	300 ng/mL*	CE
DPP-114	Urine	Panel	25 T		CE
DPP-101	Urine	Dipstick	50 T		CE

Synthetic Marijuana (K2)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DSM-102	Urine	Cassette	25 T	50 ng/mL*, 30 ng/mL, 25 ng/mL	CE
DSM-114	Urine	Panel	25 T		CE
DSM-101	Urine	Dipstick	50 T		CE

Tricyclic Antidepressants (TCA)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DTC-102	Urine	Cassette	25 T	1000 ng/mL*	CE
DTC-114	Urine	Panel	25 T		CE
DTC-101	Urine	Dipstick	50 T		CE

Marijuana (THC) FDA 510(k) Clearances

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DTH-102	Urine	Cassette	25 T	25 ng/mL, 20 ng/mL, 50 ng/mL*, 150 ng/mL, 300 ng/mL, 600 ng/mL	CE
DTH-114	Urine	Panel	25 T		CE
DTH-101	Urine	Dipstick	50 T		CE

Tramadol (TML)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DTM-102	Urine	Cassette	25 T	100 ng/mL*, 300 ng/mL	CE
DTM-114	Urine	Panel	25 T		CE
DTM-101	Urine	Dipstick	50 T		CE

Zolpidem (ZOL)

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DZO-102	Urine	Cassette	25 T	50 ng/mL*	CE
DZO-114	Urine	Panel	25 T		CE
DZO-101	Urine	Dipstick	50 T		CE

Multi-Drug 2~12 Drugs Rapid Test Cassette

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DOA-125~1125	Urine	Cassette	25 T	Ref: Analytes	CE

Multi-Drug 2~12 Drugs Rapid Test Panel

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DOA-124~1124	Urine	Panel	25 T	Ref: Analytes	CE

Multi-Drug 2~12 Drugs Rapid Test 1-Step Cup

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DOA-127~1127-A	Urine	Cup	25 T	Ref: Analytes	CE

Multi-Drug 2~12 Drugs Rapid Test 2-Step Cup

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DOA-127~1127-B	Urine	Cup	25 T	Ref: Analytes	CE

Adulteration Strip

Cat. No.	Specimen	Format	Kit Size	Analytes	CE Status
DAT-101	Urine	Dipstick	25 T	OX, NIT, SG, GLU, pH, CRE	CE

* Standard Cut-off



Drugs of Abuse

Specimen: Oral Fluid

Cat. No.	Analyte	Format	Kit Size	Cut-Off	CE Status
DAM-802	Amphetamine (AMP)	Cassette	25 T	50 ng/mL	CE
DBZ-802	Benzodiazepines (BZO)	Cassette	25 T	20 ng/mL	CE
DCO-802	Cocaine (COC)	Cassette	25 T	20 ng/mL	CE
DCT-802	Cotinine (COT)	Cassette	25 T	20 ng/mL	CE
DMA-802	6-Monoacetylmorphine (6-MAM)	Cassette	25 T	10 ng/mL	CE
DMD-802	Ecstasy (MDMA)	Cassette	25 T	50 ng/mL	CE
DME-802	Methamphetamine (MET)	Cassette	25 T	50 ng/mL	CE
DMT-802	Methadone (MTD)	Cassette	25 T	30 ng/mL	CE
DOP-802	Opiates (OPI)	Cassette	25 T	40 ng/mL	CE
DOX-802	Oxycodone (OXY)	Cassette	25 T	20 ng/mL	CE
DPC-802	Phencyclidine (PCP)	Cassette	25 T	10 ng/mL	CE
DSM-802	Synthetic Marijuana (K2)	Cassette	25 T	25 ng/mL	CE
DTH-802	Marijuana (THC)	Cassette	25 T	20 ng/mL	CE

Multi-Drug 2~6 Drugs Rapid Test Midstream

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DSD-823~863	Oral Fluid	Midstream	25 T	Ref: Analytes	CE

Multi-Drug 2~6 Drugs Rapid Test Cassette

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DSD-825~865	Oral Fluid	Cassette	25 T	Ref: Analytes	CE

Multi-Drug 2~10 Drugs Rapid Test Cup

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DSD-827~8107	Oral Fluid	Cup	25 T	Ref: Analytes	CE

Specimen: Powder (Forensic Use)

Cat. No.	Analyte	Format	Kit Size	Cut-Off	CE Status
DAC-X14	Acetaminophen (ACE)	Panel	25 T	5000 ng/mL	CE
DAM-X14	Amphetamine (AMP)	Panel	25 T	1000 ng/mL	CE
DBA-X14	Barbiturate (BAR)	Panel	25 T	300 ng/mL	CE
DBU-X14	Buprenorphine (BUP)	Panel	25 T	10 ng/mL	CE
DBZ-X14	Benzodiazepines (BZO)	Panel	25 T	300 ng/mL	CE
DCO-X14	Cocaine (COC)	Panel	25 T	300 ng/mL	CE
DCT-X14	Cotinine (COT)	Panel	25 T	200 ng/mL	CE
DDI-X14	Diazepam (DIA)	Panel	25 T	300 ng/mL	CE
DFY-X14	Fentanyl (FYL)	Panel	25 T	200 ng/mL	CE
DKE-X14	Ketamine (KET)	Panel	25 T	1000 ng/mL	CE
DMA-X14	6-Monoacetylmorphine (6-MAM)	Panel	25 T	10 ng/mL	CE
DMD-X14	Ecstasy (MDMA)	Panel	25 T	500 ng/mL	CE
DMDA-X14	Tenamfetamine (MDA)	Panel	25 T	500 ng/mL	CE
DMP-X14	Methylphenidate(MPD)	Panel	25 T	300 ng/mL	CE
DME-X14	Methamphetamine (MET)	Panel	25 T	1000 ng/mL	CE
DMO-X14	Morphine (MOP)	Panel	25 T	300 ng/mL	CE
DMQ-X14	Methaqualone (MQL)	Panel	25 T	300 ng/mL	CE
DMT-X14	Methadone (MTD)	Panel	25 T	300 ng/mL	CE
DOX-X14	Oxycodone (OXY)	Panel	25 T	100 ng/mL	CE
DPC-X14	Phencyclidine (PCP)	Panel	25 T	25 ng/mL	CE
DPP-X14	Propoxyphene (PPX)	Panel	25 T	300 ng/mL	CE
DTC-X14	Tricyclic Antidepressants (TCA)	Panel	25 T	1000 ng/mL	CE
DTM-X14	Tramadol (TML)	Panel	25 T	100 ng/mL	CE
DLS-X14	Lysergic Acid Diethylamide (LSD)	Panel	25 T	20 ng/mL	CE
DZO-X14	Zolpidem(ZOL)	Panel	25 T	50 ng/mL	CE

Multi-Drug 2~12 Drugs Rapid Test Panel

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
DOA-X24~X124	Powder	Panel	25 T	Ref: Analytes	CE

Cardiac Markers



CK-MB

Cat. No.	Specimen	Format	Kit Size	CE Status
CCK-402	WB / S / P	Cassette	25 T	CE

C-reactive Protein

Cat. No.	Specimen	Format	Kit Size	CE Status
CCR-402	WB / S / P	Cassette	20 T	CE
CCR-U401	WB / S / P	Dipstick	20 T	CE
CCR-T402	WB / S / P	Cassette	20 T	CE

Cardiac Troponin I (cTnI)

Cat. No.	Specimen	Format	Kit Size	CE Status
CTI-402	WB / S / P	Cassette	25 T	CE

D-dimer

Cat. No.	Specimen	Format	Kit Size	CE Status
CDD-402	WB / P	Cassette	25 T	CE

H-FABP

Cat. No.	Specimen	Format	Kit Size	CE Status
CFA-402	WB / S / P	Cassette	25 T	CE

H-FABP and cTnI Combo

Cat. No.	Specimen	Format	Kit Size	CE Status
CFC-425	WB / S / P	Cassette	25 T	CE

Myoglobin/CK-MB/Troponin I Combo

Cat. No.	Specimen	Format	Kit Size	CE Status
CMA-435	WB / S / P	Cassette	25 T	CE

Myoglobin

Cat. No.	Specimen	Format	Kit Size	CE Status
CMY-402	WB / S / P	Cassette	25 T	CE

Procalcitonin (PCT)

Cat. No.	Specimen	Format	Kit Size	CE Status
CPC-402	WB / S / P	Cassette	25 T	CE

5

Tumor Markers

AFP

Cat. No.	Specimen	Format	Kit Size	CE Status
TAF-402	WB / S / P	Cassette	25 T	CE

CEA

Cat. No.	Specimen	Format	Kit Size	CE Status
TCE-402	WB / S / P	Cassette	25 T	CE

FOB

Cat. No.	Specimen	Format	Kit Size	CE Status
TFO-601	Feces	Dipstick	25 T	CE
TFO-602	Feces	Cassette	25 T	CE

Hb/Hp-Hb Combo

Cat. No.	Specimen	Format	Kit Size	CE Status
THPB-625	Feces	Cassette	25 T	CE

Transferrin and FOB Combo

Cat. No.	Specimen	Format	Kit Size	CE Status
TTFC-625	Feces	Cassette	25 T	CE

PSA

Cat. No.	Specimen	Format	Kit Size	CE Status
TPS-402	WB / S / P	Cassette	25 T	Non-CE



Women's Health

hCG Pregnancy Tests (For Professional Use) FDA 510(k) Clearances

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
FHC-101	Urine	Dipstick	50 T	20 mIU/mL, 25 mIU/mL*	CE
FHC-U101	Urine	Dipstick	50 T	10 mIU/mL	CE
FHC-111+	Urine	Dipstick	25 T	20 mIU/mL, 25 mIU/mL*	CE
FHC-U111+	Urine	Dipstick	25 T	10 mIU/mL	CE
FHC-201	S / P / U	Dipstick	50 T	20 mIU/mL, 25 mIU/mL*	CE
FHC-U201	S / P / U	Dipstick	50 T	10 mIU/mL	CE
FHC-211+	S / P / U	Dipstick	25 T	20 mIU/mL, 25 mIU/mL*	CE
FHC-U211+	S / P / U	Dipstick	25 T	10 mIU/mL	CE
FHC-102	Urine	Cassette	25 T	20 mIU/mL, 25 mIU/mL*	CE
FHC-U102	Urine	Cassette	25 T	10 mIU/mL	CE
FHC-202	S / P / U	Cassette	25 T	20 mIU/mL, 25 mIU/mL*	CE
FHC-U202	S / P / U	Cassette	25 T	10 mIU/mL	CE
FHC-U402	WB	Cassette	25 T	10 mIU/mL	CE

hCG Pregnancy Tests (For Self-Testing) FDA 510(k) Clearances

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
FHC-101H	Urine	Dipstick	1 T	20 mIU/mL, 25 mIU/mL*	CE0197
FHC-U101H	Urine	Dipstick	1 T	10 mIU/mL	CE0197
FHC-102H	Urine	Cassette	1 T	20 mIU/mL, 25 mIU/mL*	CE0197
FHC-U102H	Urine	Cassette	1 T	10 mIU/mL	CE0197
FHC-103	Urine	Midstream	No Box	20 mIU/mL, 25 mIU/mL*	CE0197
FHC-103H	Urine	Midstream	1 T	20 mIU/mL, 25 mIU/mL*	CE0197
FHC-U103	Urine	Midstream	No Box	10 mIU/mL	CE0197
FHC-U103H	Urine	Midstream	1 T	10 mIU/mL	CE0197

LH Ovulation Tests (For Professional Use) FDA 510(k) Clearances

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
FLH-101	Urine	Dipstick	50 T	35 mIU/mL, 40 mIU/mL*	CE
FLH-111+	Urine	Dipstick	25 T	35 mIU/mL, 40 mIU/mL*	CE
FLH-U101	Urine	Dipstick	50 T	30 mIU/mL	CE
FLH-U111+	Urine	Dipstick	25 T	30 mIU/mL	CE
FLH-102	Urine	Cassette	25 T	35 mIU/mL, 40 mIU/mL*	CE
FLH-U102	Urine	Cassette	25 T	30 mIU/mL	CE

LH Ovulation Tests (For Self-Testing) FDA 510(k) Clearances

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
FLH-101H	Urine	Dipstick	5 T	35 mIU/mL, 40 mIU/mL*	CE0197
FLH-U101H	Urine	Dipstick	5 T	30 mIU/mL	CE0197
FLH-102H	Urine	Cassette	5 T	35 mIU/mL, 40 mIU/mL*	CE0197
FLH-U102H	Urine	Cassette	5 T	30 mIU/mL	CE0197
FLH-103	Urine	Midstream	No Box	35 mIU/mL, 40 mIU/mL*	CE0197
FLH-103H	Urine	Midstream	5 T	35 mIU/mL, 40 mIU/mL*	CE0197
FLH-U103	Urine	Midstream	No Box	30 mIU/mL	CE0197
FLH-U103H	Urine	Midstream	5 T	30 mIU/mL	CE0197

Fetal Fibronectin

Cat. No.	Specimen	Format	Kit Size	CE Status
FFF-502	Vaginal Discharge	Cassette	25T	CE

Insulin-like Growth Factor-binding Protein 1

Cat. No.	Specimen	Format	Kit Size	CE Status
FIG-502	Vaginal Discharge	Cassette	25T	CE

Toxo (IgG, IgM, IgG & IgM)

Cat. No.	Specimen	Format	Kit Size	CE Status
ITG-302	S / P	Cassette	25 T	Non-CE
ITM-302	S / P	Cassette	25 T	Non-CE
ITGM-302	S / P	Cassette	25 T	Non-CE
ITGM-325	S / P	Cassette	25 T	Non-CE

Rubella (IgG, IgM, IgG & IgM)

Cat. No.	Specimen	Format	Kit Size	CE Status
IRG-302	S / P	Cassette	25 T	Non-CE
IRM-302	S / P	Cassette	25 T	Non-CE
IRGM-325	S / P	Cassette	25 T	Non-CE

CMV (IgG, IgM, IgG & IgM)

Cat. No.	Specimen	Format	Kit Size	CE Status
ICG-302	S / P	Cassette	25 T	Non-CE
ICM-302	S / P	Cassette	25 T	Non-CE
ICGM-325	S / P	Cassette	25 T	Non-CE

+ Canister Pack

* Standard Cut-off

HSV 1/2 (IgG, IgM, IgG & IgM)

Cat. No.	Specimen	Format	Kit Size	CE Status
ISG-302	S / P	Cassette	25 T	CE
ISM-302	S / P	Cassette	25 T	CE
ISGM-325	S / P	Cassette	25 T	CE

ToRCH IgM Combo (Toxo, Rubella, CMV, HSV 1/2)

Cat. No.	Specimen	Format	Kit Size	CE Status
ITOM-345	S / P	Cassette	25 T	Non-CE

Vaginal pH

Cat. No.	Specimen	Format	Kit Size	CE Status
IVPH-504	Vaginal Discharge	Panel	25 T	CE

Chlamydia & Gonorrhea

Cat. No.	Specimen	Format	Kit Size	CE Status
ICH-502	Swab / Urine	Cassette	20 T	Non-CE
IGO-502	Swab	Cassette	20 T	Non-CE

Infectious Disease

Adenovirus & Rotavirus

Cat. No.	Specimen	Format	Kit Size	CE Status
IAD-602	Feces	Cassette	25 T	CE
IRO-602	Feces	Cassette	25 T	CE
IRAC-625	Feces	Cassette	25 T	CE

Chagas IgG

Cat. No.	Specimen	Format	Kit Size	CE Status
ICHA-402	WB / S / P	Cassette	25 T	CE

Chikungunya IgG & IgM

Cat. No.	Specimen	Format	Kit Size	CE Status
ICHI-402	WB / S / P	Cassette	25 T	CE

Dengue (IgG & IgM, NS1 Ag, Ag & Ab Combo)

Cat. No.	Specimen	Format	Kit Size	CE Status
IDEN-402	WB / S / P	Cassette	25 T	CE
IDES-402	WB / S / P	Cassette	25 T	CE
IDEC-425	WB / S / P	Cassette	25 T	CE

Filariasis IgG&IgM

Cat. No.	Specimen	Format	Kit Size	CE Status
IFIL-302	S / P	Cassette	25 T	CE

Hepatitis A IgM

Cat. No.	Specimen	Format	Kit Size	CE Status
IHA-302	S / P	Cassette	25 T	CE

Hepatitis B (HBcAb, HBsAb, HBsAg, HBsAb, HBsAg, HBV Combo)

Cat. No.	Specimen	Format	Kit Size	CE Status
IHBcb-302	S / P	Cassette	25 T	CE
IHBeb-302	S / P	Cassette	25 T	CE
IHBeg-302	S / P	Cassette	25 T	CE
IHBsb-302	S / P	Cassette	25 T	CE
IHBsg-301	S / P	Dipstick	50 T	Non-CE
IHBsg-302	S / P	Cassette	25 T	Non-CE
IHBsg-401	WB / S / P	Dipstick	50 T	Non-CE
IHBsg-402	WB / S / P	Cassette	25 T	Non-CE
IHB-355	S / P	Cassette	25 T	Non-CE

Hepatitis C

Cat. No.	Specimen	Format	Kit Size	CE Status
IHC-302	S / P	Cassette	25 T	Non-CE
IHC-402	WB / S / P	Cassette	25 T	Non-CE

Hepatitis E IgG & IgM

Cat. No.	Specimen	Format	Kit Size	CE Status
IHE-302	S / P	Cassette	25 T	CE



Infectious Disease

HIV 1.2.O

Cat. No.	Specimen	Format	Kit Size	CE Status
IHI-T402	WB / S / P	Cassette	25 T	Non-CE

H. pylori Antibody

Cat. No.	Specimen	Format	Kit Size	CE Status
IHP-302	S / P	Cassette	25 T	CE
IHP-402	WB / S / P	Cassette	25 T	CE

H. pylori Antigen

Cat. No.	Specimen	Format	Kit Size	CE Status
IHP-602	Feces	Cassette	25 T	CE

Influenza A

Cat. No.	Specimen	Format	Kit Size	CE Status
IFLA-502	Swab / Nasal Aspirate	Cassette	20 T	CE

Influenza A+B

Cat. No.	Specimen	Format	Kit Size	CE Status
IIN-502	Swab / Nasal Aspirate	Cassette	20 T	CE

Leishmania IgG & IgM

Cat. No.	Specimen	Format	Kit Size	CE Status
ILEI-402	WB / S / P	Cassette	25 T	CE

Leptospira IgG & IgM

Cat. No.	Specimen	Format	Kit Size	CE Status
ILEP-402	WB / S / P	Cassette	25 T	CE

Malaria (P.f. , P.f./Pan , P.f./P.v.)

Cat. No.	Specimen	Format	Kit Size	CE Status
IMA-402	WB	Cassette	25 T	CE
IMPN-402	WB	Cassette	25 T	CE
IMPV-402	WB	Cassette	25 T	CE

MONO

Cat. No.	Specimen	Format	Kit Size	CE Status
IMO-402	WB / S / P	Cassette	25 T	CE

Norovirus

Cat. No.	Specimen	Format	Kit Size	CE Status
INO-602	Feces	Cassette	25 T	CE

Strep A

Cat. No.	Specimen	Format	Kit Size	CE Status
IST-501	Throat Swab	Dipstick	25 T	CE
IST-502	Throat Swab	Cassette	20 T	CE
ISTB-501	Throat Swab	Dipstick	25 T	CE
ISTB-502	Throat Swab	Cassette	20 T	CE

Strep B

Cat. No.	Specimen	Format	Kit Size	CE Status
ISB-501	Swab	Dipstick	25 T	CE
ISB-502	Swab	Cassette	20 T	CE

Syphilis IgG & IgM

Cat. No.	Specimen	Format	Kit Size	CE Status
ISY-301	S / P	Dipstick	50 T	CE
ISY-302	S / P	Cassette	25 T	CE
ISY-401	WB / S / P	Dipstick	50 T	CE
ISY-402	WB / S / P	Cassette	25 T	CE

Tetanus

Cat. No.	Specimen	Format	Kit Size	CE Status
ITE-402	WB / S / P	Cassette	25 T	CE

Tuberculosis

Cat. No.	Specimen	Format	Kit Size	CE Status
ITB-402	WB / S / P	Cassette	25 T	CE

Typhoid IgG & IgM

Cat. No.	Specimen	Format	Kit Size	CE Status
ITY-402	WB / S / P	Cassette	25 T	CE

Others



Blood Stain

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
OBS-901	Stain	Dipstick	25 T	50ng/ml	CE
OBS-902	Stain	Cassette	20 T	50ng/ml	CE

Human Semen

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
OHS-901	Swab	Dipstick	25 T	4ng/ml	CE
OHS-902	Swab	Cassette	20 T	4ng/ml	CE

Micro-Albumin Semi-Quantitative

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
OMIA-101	Urine	Dipstick	50 T	See Insert	CE

Micro-Albumin Qualitative

Cat. No.	Specimen	Format	Kit Size	Cut-Off	CE Status
OMAL-101	Urine	Dipstick	50 T	20µg/ml	CE



Acro Biotech is a California-based manufacturer of rapid tests offering the products from state-of-the-art ISO13485:2003 accredited CMDCAS-compliant manufacturing facilities. Our core team comprises of seasoned professionals in the field of rapid tests with more than two decades of experience in reputed rapid test manufacturing companies, working in senior positions. Our product range includes USFDA (510)K approved DOA Tests as well as an expansive range comprising of Women's Health related rapid tests, Cardiac Markers, Tumor Markers and Infectious Disease Tests for International Market.

Our Tests work with a wide range of specimens, including urine, oral Fluid, blood, serum, plasma and vaginal swabs.

At The Cusp of Innovation and Commercialization, we are working hard to introduce new products in the market that would set the future growth path, including the combination of digitization and ease of rapid tests with our forthcoming introductions of Lateral Flow Test Readers.

9500 Seventh Street, Unit M, Rancho Cucamonga,
CA 91730, U.S.A.
Tel: (909) 466-6892

www.acrobiotech.com
info@acrobiotech.com