



GENERAL INFORMATION

Throughput	Up to 600 tests/hour for per module
Fastest time to first result	12 minutes
Detection principle	Magnetic particle-based direct chemiluminescence
Chemiluminescence label	Acridinium Ester
Test Method	Sandwich, competitive, indirect



SAMPLE LOADING SYSTEM

Sample position	140/480 sample positions with STAT function
Sample tube	Microcup and blood collecting tube
Sample volume	5-150 µL
Reaction vessels loader capacity	3000 pcs, with continuous loading



REAGENT SYSTEM

Reagent position	40 reagent positions with on-board refrigeration
Reagent carousel temperature	2-8°C
Reagent detection	Reagent level tracking
Reagent pack	50/100/150* tests per kit



REACTION SYSTEM

Capacity	230 reaction vessels
Incubation temperature	37.0°C ± 0.3°C
Mixing system	Non-contact vortex mixer



MAGNETIC SEPARATION SYSTEM

Washing method	4-step magnetic separation
Temperature	37.0 °C ± 0.3 °C



OPERATING CONDITIONS

Power supply	100-240 V, 50/60 Hz
Power consumption	1200 VA
Temperature	10-30 °C
Ambient humidity	20%-85% RH
Atmospheric pressure	61.6-106.0 kPa
Dimensions & Weight	Main analyzer (without sample loading unit) 1130 * 1050 * 1460mm (L*W*H), 390kg ASU sample loading unit: 350 * 1040 * 900mm (L*W*H), 70kg iLoader-SU sample loading unit: 880 * 1040 * 1690mm (L*W*H), 230kg

Local Distributor:

P/N: ENG-i302-202503



* Availability may vary by region. Please confirm with your YHLO sales representative.



YHLO Stock code
SSE: 688575
SHENZHEN YHLO BIOTECH CO., LTD.

SHENZHEN YHLO BIOTECH CO., LTD.
Building 1, YHLO Biopark
Baolong 2nd Road, Baolong Subdistrict,
Longgang District,
Shenzhen
518116, Guangdong,
China
Tel: +86 755 26601910
Email: marketing@szyhlo.com
Website: www.szyhlo.com

YHLO

FULLY AUTOMATED
CHEMILUMINESCENCE IMMUNOASSAY ANALYZER

iFlash
9000



iFlash 9000

Flexible

With STAT positions

Support up to 4-module integration with iFlash 9000, iBC biochemistry series and iTLA system

Continuous loading of samples and reagents



Intelligent

Reflex, auto-dilution and auto-rerun functions

Real-time monitoring of reagents and consumables

Probes featured with collision prevention, and liquid level, empty aspiration and clot detection functions via capacitive sensors




Liberate

Capacity of 3000 RV with auto-replenishment system

Reagent carousel with automatic micro-magnetic bead mixing function

iFlash 9000



Stable

Array-type surrounding magnetic field design, precisely controls washing procedure for optimal repeatability

Non-contact vortex mixer, minimizing carry-over risk

Automatic routine maintenance



Accurate

Temperature-controlled, incubation, magnetic separation and reagent carousels, fluctuations within $\pm 0.1^{\circ}\text{C}$

Probes with high-pressure inner & outer walls washing and precise aspiration with an accuracy of $\pm 0.5\ \mu\text{L}$



Harmonious

Compact and user-friendly design adapt to diverse laboratory environments

Offers both concentrated and lyophilized* wash buffers for easy operation

PRODUCT SPECIFICATIONS

Throughput:

UP TO 600 TESTS 1 HOUR / MODULE

Optional sample positions:

ASU: 140 ILOADER: 480

Integration:

UP TO 4 MODULES



High-Throughput CLIA Analyzer

- Simultaneous processing of large-scale test volume
- Fastest first result time: 12 min
- High space efficiency: 386 T/H/m2



Designed for Large Central Labs

- Large onboard capacity for consumables, extended walk away time
- Supports integration with iBC biochemistry series or iTLA total laboratory automation system

iFlash 9000

Chemiluminescence Immunoassay Analyzer

Routine		Reproductive Health	Autoimmunity	Infection
Tumor Marker PG I PG II AFP CEA CA 125 CA 15-3 CA 19-9 Free PSA Total PSA CYFRA 21-1 NSE SCC HE4 CA72-4 Gastrin-17 ProGRP CA242 CA50**	Growth Hormone hGH IGF-1** IGFBP-3**	Fertility AMH Inhibin B Inhibin A LH FSH HCG Prolactin DHEA-S Unconjugated Estriol E2 PAPP-A free βhCG Progesterone Testosterone 17α-OH-P Androstenedione SHBG	Rheumatoid Arthritis Anti-CCP RF IgG RF IgM RF RA33 IgG	Blood Bank Anti-TP (Syphilis) Anti-HCV HBsAg Anti-HBs HBeAg Anti-HBe Anti-HBc HIV Combo
	Diabetes Insulin C-Peptide Adiponectin Glucagon**		Connective Tissue Disease ANA dsDNA IgG Sm IgG SS-A IgG SS-B IgG Scl-70 IgG Jo-1 IgG RNP70 IgG Rib-P0 IgG CenpB IgG PM-Scl IgG	
	Cardiac Marker CK-MB Myoglobin Troponin-I hs-cTnI BNP NT-proBNP D-Dimer		Autoimmune Liver SMA IgG AMA-M2	
Thyroid T4 FT4 T3 FT3 TSH Tg Anti-Tg Anti-TPO Anti-TSHR Calcitonin**	Hypertension ACTH Cortisol Renin** Aldosterone**	TORCH Toxo IgG Toxo IgM CMV IgG CMV IgM HSV-1 IgG HSV-2 IgG HSV-1 IgM HSV-2 IgM Rubella IgG Rubella IgM	Vasculitis Diagnostics PR3 IgG GBM IgG MPO IgG	Respiratory Disease Mycoplasma pn. IgG Mycoplasma pn. IgM Chlamydia pn. IgG Chlamydia pn. IgM
	Bone Metabolism 25-OH Vitamin D PTH Calcitonin** Osteocalcin**		Autoimmune Diabetes IA-2A IAA GADA ICA ZnT8A	
Thrombosis Marker TAT** PIC** TM** t-PAIC** D-Dimer FDP**	Liver Fibrosis HA PIIIPN-P Col IV Laminin	Anti-Phospholipid Syndrome Cardiolipin IgG Cardiolipin IgM Cardiolipin IgA Anti-Cardiolipin β2-Glycoprotein I IgG β2-Glycoprotein I IgA β2-Glycoprotein I IgM Anti-β2-Glycoprotein I	Celiac Disease tTg IgA tTg IgG DGP IgA DGP IgG	EBV EB VCA IgG EB VCA IgM EB VCA IgA EB EA IgM EB NA IgG EB NA IgA
Anemia Ferritin Folate Vitamin B12			Allergy Total IgE	Inflammation PCT IL-6 CRP**
			IgG Subclass IgG1 IgG2 IgG3 IgG4	

** Coming Soon